


IN-BIOS[9,5]000 2020

Illumina Technology

Arvind Sundaram
Oct 19, 2020

Norwegian Sequencing Centre
OUS, Ullevål, Oslo



DNA sequencing

- ❖ First generation - **past, present**
 - ❖ Up to 1 kb; high quality data; multiplexed
 - ❖ SANGER; Highly automated (ABI Sanger 3730xl)
- ❖ Second generation - **present**
 - ❖ Shorter reads; Massive parallelisation and real high throughput
 - ❖ Illumina, BGISEq, Ion-torrent, [454, Solid]
 - ❖ RNA is reverse-transcribed to cDNA before sequencing
- ❖ Third generation - **[present] future**
 - ❖ Long-read sequencing; Single-molecule sequencing (without amplification)
 - ❖ PacBio, Oxford Nanopore, [more in development]
 - ❖ Potential to sequence RNA directly

High throughput sequencing

illumina[®]

MiniSeq

MiSeq

NextSeq

HiSeq 2500

HiSeq 3 / 4000

HiSeq X

NovoSeq

BGI **SEQ**

Roche 454

SOLiD

Ion Torrent



PACBIO[®]

RS II

Sequel



MinION

Flongle

GridION

PromethION

Illumina sequencers



Benchtop sequencers

Production-scale sequencers


Data output: 144 Mb - 500 Gb
Read length: 25 - 300 nt
Read type: Single / Paired end

9 Gb - 2400 Gb
50 - 250 nt
Single / Paired end




Illumina sequencers

- ❖ Second generation sequencing technique
- ❖ Sequencing-by-synthesis aka SBS
 - ❖ <https://www.youtube.com/watch?v=fCd6B5HRaZ8>
- ❖ Mass parallelisation and real high throughput

Illumina sequencers




Sequencing



Library prep and sequencing


Fragment (DNA) sequenced: up to 800 bp

A. Library Preparation




NGS library is prepared by fragmenting a gDNA sample and ligating specialized adapters to both fragment ends.

B. Cluster Amplification



Library is loaded into a flow cell and the fragments are hybridized to the flow cell surface. Each bound fragment is amplified into a clonal cluster through bridge amplification.

C. Sequencing



Sequencing reagents, including fluorescently labeled nucleotides, are added and the first base is incorporated. The flow cell is imaged and the emission from each cluster is recorded. The emission wavelength and intensity are used to identify the base. This cycle is repeated "n" times to create a read length of "n" bases.

Library prep and sequencing

Fragment (DNA) sequenced: up to 800 bp



Add adapters during library preparation


Multiplexing: single / dual index



Multiplexing: pooling



Read type: Single end



Read type: Paired end



Library prep

<https://www.illumina.com/techniques/sequencing/ngs-library-prep/library-prep-methods.html>

Try the “Sequencing
Methods selector”

Check out 3 posters: DNA,
RNA and single cell



Single/dual indexed samples

MiSeq
HiSeq 2500
NovaSeq

Figure 1 Single-Indexed Sequencing




Figure 2 Dual-Indexed Sequencing on a Paired-End Flow Cell (Workflow A)





Figure 5 Dual-Indexed Sequencing on a Single-Read Flow Cell (HiSeq 2500 or HiSeq 2000)



Single/dual indexed samples

iSeq
MiniSeq
NextSeq
HiSeq 3/4000
HiSeq X

Figure 1 Single-Indexed Sequencing






Figure 3 Dual-Indexed Sequencing on a Paired-End Flow Cell (Workflow B)



Figure 4 Dual-Indexed Sequencing on a Single-Read Flow Cell (HiSeq 4000 or HiSeq 3000)





Four vs two colour chemistry



Old machines use 4 colours
New machines use 2 colours

Sequencing is almost
twice as fast



Illumina sequencers



Benchtop sequencers

Production-scale sequencers

Read length:
Single end

Read length:
Paired end

Data output:

36

25

100

250

144 Mb

50

50

125

300

-


75

75


150

2400 Gb


Data output



Data output



Known issues



What can you sequence using Illumina

- ❖ DNA studies
 - ❖ Whole genome sequencing - short reads are a pitfall
 - ❖ Genome re-sequencing
 - ❖ Exomes and target re-sequencing...
 - ❖ ChiP seq and more...
- ❖ RNA studies
- ❖ modification studies
 - ❖ Methylation and more...



<https://www.illumina.com/techniques/sequencing/ngs-library-prep/library-prep-methods.html>