

*INF-BIOx121 2017*

---

HTS

Arvind Sundaram

Sep 04, 2017

Norwegian Sequencing Centre  
OUS, Ullevål, Oslo


---



# High throughput sequencing

[illegible]

# High throughput sequencing



---

# High throughput (DNA) sequencing




---

- ❖ First generation
- ❖ Second generation
- ❖ Third generation
- ❖ What about RNA???
- ❖ unstable
- ❖ amplification issues
- ❖ reverse-transcribed to cDNA before being sequenced



# Sequencing

Dideoxynucleotides  
ddNTPs



---

# First generation DNA sequencing

---

- ❖ SANGER
- ❖ Highly automated (ABI Sanger 3730xl)
- ❖ Up to 1 kb; high quality data; multiplexed

---

# Second generation DNA sequencing

---

- ❖ Sequencing-by-synthesis aka SBS
  - ❖ <https://www.youtube.com/watch?v=fCd6B5HRaZ8>
- ❖ Mass parallelisation and real high throughput

---

# Third generation DNA sequencing

---

- ❖ Single molecule sequencing
- ❖ Real time sequencing
- ❖ Non-SBS methods
- ❖ Mostly under-development
- ❖ Few commercially available / successful methods
  - ❖ PacBio (<https://www.youtube.com/watch?v=NHCJ8PtYCFc>)
  - ❖ Oxford Nanopore (<https://www.youtube.com/watch?v=hs0FdiTHMbc>)



# High throughput sequencing

illumina®

MiniSeq

MiSeq

NextSeq

HiSeq 2500

HiSeq 3000 / 4000

HiSeq X

NovoSeq

Roche 454

SOLiD

HeliScope



RS II

Sequel



PGM

Proton

S5 / S5XL

BGISEQ


---

# Special types


---

- ❖ 10X genomics
- ❖ Dovetail genomics
- ❖ BioNano genomics
- ❖ Moleculo / TruSeq synthetic reads

# High throughput sequencing




# High throughput sequencing





# High throughput sequencing





---

# Experimental design

---

- ❖ Platform choice
- ❖ Technology variation
  - ❖ Technical bias
    - ❖ Lane bias
    - ❖ Index / barcode bias
    - ❖ Batch effect
    - ❖ Duplicates
- ❖ Biological question
- ❖ Sample variation
- ❖ Sequencing depth
- ❖ Data analysis
- ❖ Species-specific information
  - ❖ Is there a genome sequence available??
  - ❖ Genome size (c-value) (<http://genomesize.com/>)


---

# Experimental design


---

- ❖ Short or long fragments
- ❖ Single or paired end
- ❖ Depth required
- ❖ Coverage required
- ❖ Multiplexing
  - ❖ single index
  - ❖ dual index
  - ❖ or more??
- ❖ Replicates
  - ❖ Biological
  - ❖ Technical

# High throughput sequencing





# Sequencing depth and coverage




GAllx, Genome Analyzer Iix; PacBio, Pacific Biosciences; PGM, personal genome machine.

# Sequencing depth and coverage





# Sequencing depth and coverage



Techniques	Read counts in representative studies
DNaseI-seq and FAIRE-seq	20–50 million
CLIP-seq	7.5 million; 36 million
iCLIP and PAR-CLIP	8 million; 14 million
ChIRP and CHART	26 million
4C	1–2 million
ChIA-PET	20 million
5C	25 million
Hi-C	>100 million
MeDIP-seq	60 million
CAP-seq	>20 million
ChIP-seq	>10 million per sample (point source); >20 million per sample (broad source)